

# Retract-O-Tape®

## Air-Cushioned Vascular Loops

Effective occlusion of arteries and vessels

Occludes or retracts a vessel without tissue damage worry  
Maintains round profile while providing necessary flexibility  
Sealed silicone tubing provides an air-cushion  
Needle version offers superior control for safe, quick positioning

### Materials:

Needle: Surgical stainless steel

Tubing: Radiopaque Silicone

Latex free

Available sterile and non-sterile from Quest Medical, Inc. and from all procedure kit suppliers

1041



1071



1015



1012 | 1032



1001 | 1002



### REF

1041 Blunt needle vascular loop, specially designed, pre-attached, curved blunt tip 18G needle; thin wall  
Length: 12 in. (30.48 cm) | Diameter: 0.049 in. (1.24 mm)

1015 Retract-O-Tape vascular loop; thin wall  
Length: 18 in. (45.72 cm) | Diameter: 0.131 in. (3.32 mm)

1012 Retract-O-Tape vascular loop; thin wall  
Length: 18 in. (45.72 cm) | Diameter: 0.066 in. (1.67 mm)


1032 Retract-O-Tape vascular loop; thick wall  
Length: 18 in. (45.72 cm) | Diameter: 0.066 in. (1.67 mm)

1001 Retract-O-Tape vascular loop; thin wall  
Length: 12 in. (30.48 cm) | Diameter: 0.049 in. (1.24 mm)

1002 Retract-O-Tape vascular loop; thin wall  
Length: 12 in. (30.48 cm) | Diameter: 0.066 in. (1.67 mm)

1071 Semi-pointed needle vascular loop, specially designed, pre-attached, curved blunt tip 18G needle; thin wall  
Length: 12 in. (30.48 cm) | Diameter: 0.049 in. (1.24 mm)

*Sold in 10 packets of 2 units (Total: 20 units/case)  
except \*REF 1015, 10 packets of 1 unit (Total: 10/units/case)*



**QUEST Medical, Inc.** [www.questmedical.com](http://www.questmedical.com)

© 2014 Quest Medical, Inc. Retract-O-Tape is a trademark of Quest Medical, Inc. 2010-6

[custserv@questmedical.com](mailto:custserv@questmedical.com)

**R<sub>x</sub> ONLY**

+800.627.0226 (option 1) | +972.390.9800 (v) | +972.390.7173 (f)